

# JalSakshi Synthetic Demo Reference Card v1

## DEMO ONLY - NOT WATER ANALYSIS

Protocol SYN-COLOR-001 | Demo lot DEMO-0001 | Expiry: N/A

**SYN-A**

#3B82F6

**SYN-B**

#10B981

**SYN-C**

#F59E0B

**SYN-D**

#EF4444

### Synthetic capture procedure

1. Print at 100% on matte white paper; do not color-correct.
2. Use even indoor light, no flash, and keep all four corner marks visible.
3. Place one controlled SYN-A/B/C/D tile in the center below the swatches.
4. Capture, save offline, reopen after app restart, then sync when online.
5. A supervisor may accept/refer only the synthetic demo record.


### Required evidence log

Phone model / Android version / RAM / available storage  
Printer / paper / lighting / capture distance / timestamp  
App build hash / fixture ID / expected demo label / observed label  
Offline restart result / sync result / synthetic supervisor decision  
Never record concentration, safety, potability, diagnosis, or remedy.

Printed colors vary by printer, paper, light, and camera. Mismatch means DEMO INVALID, not a water result.

# Controlled Synthetic Sample Tiles


~~CUT OR DISPLAY FOR DEMO ONLY - NOT SCIENTIFIC REFERENCE MATERIAL~~



SYN-A

**EXPECTED DEMO LABEL: SYN-A | #3B82F6**


No concentration or potability meaning



SYN-B

**EXPECTED DEMO LABEL: SYN-B | #10B981**


No concentration or potability meaning



SYN-C

**EXPECTED DEMO LABEL: SYN-C | #F59E0B**

No concentration or potability meaning



SYN-D

**EXPECTED DEMO LABEL: SYN-D | #EF4444**

No concentration or potability meaning